

Suspected COVID-19 cases

Sampling and Packaging

Diagnostic samples for suspected cases



- 1. Upper respiratory tract sample options:
 - individual nose and throat swabs in separate collection tubes OR
 - combined nose and throat swab in one collection tube containing universal transport medium OR
 - single swab used for throat then nose OR
 - nasopharyngeal aspirate.



2. Lower respiratory tract sample in universal container (sputum) if obtainable.

If the patient is admitted, take a sample for acute serology: 5mL in either serum tube or plain (no additive) tube. For children <12 years, 1mL is acceptable.

Packaging



Ensure lids are tight and decontaminate outside



Wrap sample with absorbent material



Individually place into zip-lock bag



Seal zip-lock bag & change gloves



Add all wrapped samples into second bag



Wrap bagged samples in ample packaging & place into bio-bottle



Firmly attach lid



Ensure E28 request forms for each individual sample are OUTSIDE the bio-bottle



Place E28 request forms & bio-bottle into transport box & add security seal



Label box with 'Priority 10', you can download the labelling template at bit.ly/2vrRNxT





Courier all samples to PHE via Category B UN3373.
PHE is open to receive samples 24/7

Download the request form for COVID-19 (formerly novel coronavirus 2019-nCoV)

www.gov.uk/government/publications/testing-for-wuhan-novel-cov-2019-ncov

